

Defibrillation.
Made easy. For everyone.



Why do we need public access defibrillation?

Sudden cardiac arrest (SCA) can strike any person, anywhere, anytime. When someone collapses from SCA, immediate defibrillation and cardiopulmonary resuscitation (CPR) are essential for any chance of recovery. Immediate treatment is critical for SCA victims since survival chances decrease by about 10% with every minute without defibrillation.

Authorities agree that defibrillation should be provided as early as possible, preferably within first 3 to 5 minutes. Even with the fastest emergency medical service system, the professionals may not be able to reach the victim on time. This results in the necessity of the lifesaving defibrillation to be delivered directly on the site of emergency within first minutes and by any responder.

SURVIVAL RATE INCREASES BY

90%

IF DEFIBRILLATION OCCURS WITHIN
1 MINUTE OF VICTIM COLLAPSING

THERE ARE APPROXIMATELY

60,000

OUT OF HOSPITAL CARDIAC
ARRESTS OCCURING EVERY YEAR

IN THE UK APPROXIMATELY

270

CHILDREN DIE AFTER SUFFERING
A CARDIAC ARREST AT SCHOOL



CardiAid Semi-automatic:

One-button operation to provide easier control for the rescuer.



CardiAid Fully-automatic:

No need to press the button to deliver the shock.



What is a CardiAid AED?

The CardiAid AED (Automated External Defibrillator) was designed to be the simplest and clearest possible defibrillator for public use, making it possible to administer an electric shock with a minimum of delay to anyone suffering acute heart failure.

The CardiAid AED has been designed for public use of defibrillators in mind at places such as children's day care centres, schools, sports complexes, doctors' surgeries, veterinary and dental surgeries, or any type of premises with large work populations.

The CardiAid AED is completely safe to use. Our German manufactured defibrillator only delivers a life-saving electric shock when the analysis has determined it necessary. With Fully and Semi-Automatic CardiAid AEDs you can be safe in the knowledge that a shock is only delivered when it detects a shockable rhythm. With such features as pace maker detection, you can rely on this entirely. A loud voice giving clear instructions guide you through the resuscitation process, along with clear pictograms. A metronome provides the heart massage rate. The AED produces strong visual and auditory guidance to help support you in the deployment of our defibrillator.

Features of the CardiAid AED

CardiAid automatically switches on by simply opening the cover



Pre-connected electrodes help save time and make the process simpler for the user



“Status indicator lights” provide easy interpretation of the self-test results




CardiAid guides the user step by step with a combination of clear voice instructions, illustrations and lights for the maximum efficiency. All steps are ensured to be applied accurately even if the responder has limited knowledge or experience



CardiAid Semi Automatic works with one-button only – the shock button. The Fully Automatic works independently and has no shock button. There is no on-off button, therefore no time loss!



Suitable for use in any environment



Training & Servicing. Complete peace of mind



CardiAid Training Programme

It is important to ensure calmness and confidence when using a CardiAid AED. This is why our training offer instils confidence and builds knowledge of how to react in the event of a Sudden Cardiac Arrest. Our training includes up to 8 people at your premises and can last between 2-3 hrs.

This is available to purchase separately or can be included in our rental packages.

The course covers and includes:

- What is a Cardiac Arrest
- Basic CPR awareness training
- CardiAid AED Orientation
- Process of dealing with a suspected Cardiac Arrest
- Identifying when to give CPR
- End Of Course Test
- Certificated upon course completion

CardiAid Servicing Programme

CardiAid is a high-quality product which is designed and manufactured with the highest medical standards in mind. At CardiAid we always want to ensure your defibrillator is ready to use and remains in the very best condition. After all, you never know when you might really need it.

Regular checks and periodic servicing are of the utmost importance for any defibrillator. At CardiAid we offer a fixed cost service contract which guarantees that your CardiAid defibrillator is always in good working order.

Rental or purchase? The choice is yours.

Having a CardiAid AED can be the difference between life and death. We offer tailored solutions to meet your needs from under £1 per day right through to our complete business packages that enable you to remotely monitor your CardiAid AED.

Peace of mind rental. From under £1 per day.



Do you recognise the importance of your premises or business premises being fitted with an AED, but are put off by the purchase price? If so renting the AED is the answer.

Our rental packages make it possible for every organisation to ensure a safe environment by installing an AED. For a fixed monthly fee you can have the most advanced CardiAid AED to hand for the duration of 72 months. Pediatric servicing and replacement of the batteries and electrodes is included in the subscription fee, so you will never be faced with unexpected expenses.

Purchasing with our Price Guarantee.

We will not be beaten on price on our CardiAid bundles.



CardiLink. Connected AED safety solutions.



What is CardiLink and how does it work?

CardiLink is a system devised for organisations with multiple CardiAid AED units allowing each device to be remotely monitored. Employers can save their staff time and energy with our product monitoring. CardiLink software displays the status of your devices in a dashboard in real time and analyzes your data for the long-term monitoring and quality reports. Alarms and any occurring issues will be reported to you before the product will be used with a patient. CardiLink can also provide Data Analytics of your product for clinical studies.



Always Ready

Checks and documents
AED status



Automated QMS Compliance

Certified and online
documentation



Emergency Notifications

Sends a notification when
moved and opened



Uptime Guarantee

Via pro-active servicing
and SLA



GPS Geofence

Theft detection
movement tracking



Fleet Management

Managed solutions,
API intergration



CardiAid United Kingdom

Unit 1, Morelands Trading Estate,
Bristol Road, Gloucester, GL1 5RZ

Telephone +44 (0) 1452 937656

Email info@cardiaid.co.uk

Website www.cardiaid-aed.co.uk

CardiAid[®]
Public Access Defibrillator